

Product Datasheet

Mitogen-activated protein kinase kinase kinase antibody (Biotin) (orb241041)

Description	Mitogen-activated protein kinase kinase kinase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Mitogen-activated protein kinase kinase kinase 1 protein (1216-1493AA)
Target	Map3k1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53349
Expiration Date	12 months from date of receipt.