

Product Datasheet

Multidrug resistance protein 3 antibody (HRP) (orb241051)

Description	Multidrug resistance protein 3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Phosphatidylcholine translocator ABCB4 protein (352-708AA)
Target	Streptavidin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21440
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.