


Product Datasheet

Nidogen-1 antibody (orb241056)

Description	Nidogen-1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Nidogen-1 protein (428-665AA)
Target	Nid1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10493
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat lung tissue...