



Product Datasheet

Protein S100A9 antibody (orb241072)

Description	Protein S100A9 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Mouse Protein S100-A9 protein (2-113AA)
Target	S100a9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31725
Dilution Range	WB:1:500-1:5000, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation
analysis of Mouse
sp...



Western blot
analysis of Mouse
brain ti...