

Product Datasheet

Scavenger R cysteine-rich type 1 protein M130 antibody (orb241084)

Description	Scavenger R cysteine-rich type 1 protein M130 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Scavenger receptor cysteine-rich type 1 protein M130 protein (86-365AA)
Target	Cd163
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2VLH6
Expiration Date	12 months from date of receipt.