

Product Datasheet

Thrombospondin 2 antibody (Biotin) (orb241105)

Description	Thrombospondin 2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Thrombospondin-2 protein (19-232AA)
Target	Thbs2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q03350
Expiration Date	12 months from date of receipt.