


Product Datasheet

Translationally-controlled tumor antibody (orb241116)

Description	Translationally-controlled tumor antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Mouse Translationally-controlled tumor protein (1-172AA)
Target	Tpt1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P63028
Dilution Range	IHC-P:1:10-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of human
pl...