

Product Datasheet

Translationally-controlled tumor antibody (Biotin) (orb241117)

Description	Translationally-controlled tumor antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Translationally-controlled tumor protein (1-172AA)
Target	Tpt1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P63028
Expiration Date	12 months from date of receipt.