

Product Datasheet

Tryptase beta-2 antibody (HRP) (orb241127)

Description	Tryptase beta-2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Tryptase beta-2 protein (32-276AA)
Target	Tpsb2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21845
Expiration Date	12 months from date of receipt.