

## Product Datasheet

# WD repeat-containing and planar cell polarity effector protein fritz homolog antibody (FITC) (orb241142)

|                        |                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | WD repeat-containing and planar cell polarity effector protein fritz homolog antibody (FITC)                     |
| <b>Species/Host</b>    | Rabbit                                                                                                           |
| <b>Reactivity</b>      | Mouse                                                                                                            |
| <b>Conjugation</b>     | FITC                                                                                                             |
| <b>Immunogen</b>       | Recombinant Mouse WD repeat-containing and planar cell polarity effector protein fritz homolog protein (1-722AA) |
| <b>Target</b>          | Wdpcp                                                                                                            |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                                                   |
| <b>Form/Appearance</b> | Liquid                                                                                                           |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                    |
| <b>Note</b>            | For research use only                                                                                            |
| <b>Isotype</b>         | IgG                                                                                                              |
| <b>Clonality</b>       | Polyclonal                                                                                                       |
| <b>Uniprot ID</b>      | <b>Q8C456</b>                                                                                                    |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                  |