


Product Datasheet

C5a anaphylatoxin antibody (orb241148)

Description	C5a anaphylatoxin antibody
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Pig Complement C5a anaphylatoxin protein (1-74AA)
Target	C5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01032
Dilution Range	IHC-P:1:10-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...