


Product Datasheet

Rabies virus Matrix protein antibody (orb241168)

Description	Rabies virus Matrix protein antibody
Species/Host	Rabbit
Reactivity	Rabies virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rabies virus Matrix protein (1-202AA)
Target	M
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15200
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant Mat...