



Product Datasheet

Avpr1a protein antibody (orb241188)

Description	Avpr1a protein antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Rat Vasopressin V1a receptor protein (7-52AA)
Target	Avpr1a
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN3, PH7.2
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30560
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Mouse spleen t...