

Product Datasheet

Carbonic anhydrase 1 antibody (HRP) (orb241206)

Description	Carbonic anhydrase 1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Rat Carbonic anhydrase 1 protein (2-261AA)
Target	Ca1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	BOBNN3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.