

Product Datasheet

Cathepsin K antibody (FITC) (orb241210)

Description	Cathepsin K antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Recombinant Rat Cathepsin K protein (115-329AA)
Target	Ctsk
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	035186
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.