


Product Datasheet

Cellular tumor p53 antibody (orb241212)

Description	Cellular tumor p53 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rat Cellular tumor antigen p53 protein (1-391AA)
Target	Tp53
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10361
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...