

Product Datasheet

COL1A1 antibody (Biotin) (orb241217)

Description	COL1A1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA)
Target	Col1a1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02454
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.