

Product Datasheet

Glutathione S-transferase P antibody (HRP) (orb241239)

Description	Glutathione S-transferase P antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Rat Glutathione S-transferase P protein (2-210AA)
Target	Gstp1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04906
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.