


Product Datasheet

Guanidinoacetate N-methyltransferase antibody (orb241244)

Description	Guanidinoacetate N-methyltransferase antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rat Guanidinoacetate N-methyltransferase protein (1-236AA)
Target	Gamt
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10868
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse thymus ti...