

Product Datasheet

Integrin beta-1 antibody (Biotin) (orb241257)

Description	Integrin beta-1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Integrin beta-1 protein (21-729AA)
Target	Itgb1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49134
Expiration Date	12 months from date of receipt.