

Product Datasheet

Microtubule-associated protein tau antibody (Biotin) (orb241273)

Description	Microtubule-associated protein tau antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Microtubule-associated protein tau protein (2-752AA)
Target	Mapt
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19332
Expiration Date	12 months from date of receipt.