

Product Datasheet

Muscarinic acetylcholine R M1 antibody (Biotin) (orb241277)

Description	Muscarinic acetylcholine R M1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Muscarinic acetylcholine receptor M1 protein (210-366AA)
Target	Chrm1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08482
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.