

Product Datasheet

Myosin-binding protein C, cardiac-type antibody (orb241280)

Description	Myosin-binding protein C, cardiac-type antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rat Myosin-binding protein C, cardiac-type protein (645-864AA)
Target	Mybpc3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56741
Expiration Date	12 months from date of receipt.