

## Product Datasheet

### MAB TO HCV CORE ANTIGEN (orb24130)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | MAB TO HCV CORE ANTIGEN                            |
| <b>Species/Host</b>        | Culture                                            |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | EIA                                                |
| <b>Target</b>              | Hepatitis C Virus (HCV) Core Antigen (a.a. 80-120) |
| <b>Preservatives</b>       | NaN3                                               |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Monoclonal                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |