

Product Datasheet

Tryptase beta-2 antibody (FITC) (orb241302)

Description	Tryptase beta-2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Recombinant Rat Tryptase beta-2 protein (30-274AA)
Target	Tpsb2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50343
Expiration Date	12 months from date of receipt.