

Product Datasheet

Urokinase plasminogen activator surface R antibody (HRP) (orb241311)

Description	Urokinase plasminogen activator surface R antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Rat Urokinase plasminogen activator surface receptor ptotein (25-299AA)
Target	Plaur
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49616
Expiration Date	12 months from date of receipt.