


Product Datasheet

Vasopressin V1b R, antibody (orb241315)

Description	Vasopressin V1b R, antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rat Vasopressin V1b receptor protein (343-425AA)
Target	Avpr1b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48974
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...