

Product Datasheet

Vesicle-fusing ATPase antibody (Biotin) (orb241320)

Description	Vesicle-fusing ATPase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Vesicle-fusing ATPase protein (7-209AA)
Target	Nsf
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9QUL6
Expiration Date	12 months from date of receipt.