

Product Datasheet

Zinc transporter ZIP13 antibody (orb241323)

Description	Zinc transporter ZIP13 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rat Zinc transporter ZIP13 protein (130-233AA)
Target	Slc39a13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2M1K6
Expiration Date	12 months from date of receipt.