

Product Datasheet

IL12 p40 antibody (HRP) (orb241353)

Description	IL12 p40 antibody (HRP)
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Sheep Interleukin-12 subunit beta protein (23-327AA)
Target	IL12B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68220
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.