

Product Datasheet

Staphylococcus aureus Transposase for transposon Tn552 antibody (orb241358)

Description	Staphylococcus aureus Transposase for transposon Tn552 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus Transposase for transposon Tn552 protein (1-480AA)
Target	Transposase for transposon Tn552
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18416
Expiration Date	12 months from date of receipt.