

Product Datasheet

Streptomyces hygroscopicus Phosphinothricin N-acetyltransferase antibody (HRP) (orb241369)

Description	Streptomyces hygroscopicus Phosphinothricin N-acetyltransferase antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Streptomyces hygroscopicus Phosphinothricin N-acetyltransferase protein (1-183AA)
Target	bar
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16426
Expiration Date	12 months from date of receipt.