

Product Datasheet

Thermus aquaticus DNA polymerase I, thermostable antibody (orb241370)

Description	Thermus aquaticus DNA polymerase I, thermostable antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Thermus aquaticus DNA polymerase I, thermostable protein (310-832AA)
Target	polA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19821
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...