

Product Datasheet

Torpedo californica Acetylcholine R subunit alpha antibody (Biotin) (orb241375)

Description	Torpedo californica Acetylcholine R subunit alpha antibody (Biotin)
Species/Host	Rabbit
Reactivity	Other
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Torpedo californica Acetylcholine receptor subunit alpha protein (25-234AA)
Target	CHRNA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02710
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.