


Product Datasheet

dcyD antibody (orb241405)

Description	dcyD antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli D-cysteine desulfhydrase protein (2-328AA)
Target	dcyD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P76316
Expiration Date	12 months from date of receipt.



Western blot analysis of DH5a whole cell...