


Product Datasheet

Sus scrofa Interleukin-33 antibody (orb241417)

Description	Sus scrofa Interleukin-33 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Pig Interleukin-33 protein (123-276AA)
Target	IL-33
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	M5B263
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...