

Product Datasheet

Thrombospondin-1 antibody (FITC) (orb241419)

Description	Thrombospondin-1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Thrombospondin-1 protein (19-350AA)
Target	Thbs1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35441
Expiration Date	12 months from date of receipt.