

Product Datasheet

Crotalus adamanteus Thrombin-like enzyme crotalase antibody (orb241429)

Description	Crotalus adamanteus Thrombin-like enzyme crotalase antibody
Species/Host	Rabbit
Reactivity	Reptile
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Crotalus adamanteus Thrombin-like enzyme crotalase protein (25-262AA)
Target	Thrombin-like enzyme crotalase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	F8S114
Expiration Date	12 months from date of receipt.