

Product Datasheet

Vitellogenin-2 antibody (orb241447)

Description	Vitellogenin-2 antibody
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Gallus gallus Vitellogenin-2 protein (26-653AA)
Target	VTG2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02845
Expiration Date	12 months from date of receipt.