


Product Datasheet

Helicobacter pylori Urease subunit alpha antibody (orb241451)

Description	Helicobacter pylori Urease subunit alpha antibody
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Helicobacter pylori Urease subunit alpha protein (1-238AA)
Target	ureA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14916
Expiration Date	12 months from date of receipt.



Western blot analysis of Helicobacter py...