

Product Datasheet

52 kDa repressor of the inhibitor of the protein kinase antibody (Biotin) (orb241460)

Description	52 kDa repressor of the inhibitor of the protein kinase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human 52 kDa repressor of the inhibitor of the protein kinase protein (612-761AA)
Target	THAP12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43422
Expiration Date	12 months from date of receipt.