



Product Datasheet

ACTA1 antibody (orb241463)

Description	ACTA1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Actin, alpha skeletal muscle protein (3-377AA)
Target	Scytovirin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68133
Dilution Range	WB:1:500-5000, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Rats heart(Lane...



Western blot analysis of Mouse heart(Lan...



Immunohistochemical staining of human sk...