



Product Datasheet

Alkaline phosphatase, tissue-nonspecific isozyme antibody (orb241483)

Description	Alkaline phosphatase, tissue-nonspecific isozyme antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Alkaline phosphatase, tissue-nonspecific isozyme protein (28-501AA)
Target	ALPL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05186
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Immunohistochemical staining of human ov...