

Product Datasheet

Annexin A9 antibody (FITC) (orb241509)

Description	Annexin A9 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Annexin A9 protein (8-345AA)
Target	Defensin-like protein 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	076027
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.