

Product Datasheet

Collagen alpha-1(III) chain antibody (Biotin) (orb241564)

Description	Collagen alpha-1(III) chain antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Collagen alpha-1(III) chain protein (1222-1465AA)
Target	COL3A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02461
Expiration Date	12 months from date of receipt.