

Product Datasheet

Dual specificity protein phosphatase 13 isoform A antibody (Biotin) (orb241588)

Description	Dual specificity protein phosphatase 13 isoform A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Dual specificity protein phosphatase 13 isoform A protein (1-198AA)
Target	DUSP13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6B8I1
Expiration Date	12 months from date of receipt.