



Product Datasheet

Epidermal growth factor R antibody (orb241603)

Description	Epidermal growth factor R antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Epidermal growth factor receptor protein (26-645AA)
Target	EGFR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00533
Dilution Range	IHC-P:1:20-1:200, IF:1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunofluorescence analysis of Hela cells.