

Product Datasheet

Glucosidase 2 subunit beta antibody (FITC) (orb241633)

Description	Glucosidase 2 subunit beta antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Glucosidase 2 subunit beta protein (15-302AA)
Target	PRKCSH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14314
Expiration Date	12 months from date of receipt.