

Product Datasheet

Integrin-linked kinase-associated threonine phosphatase 2C antibody (Biotin) (orb241664)

Description	Integrin-linked kinase-associated threonine phosphatase 2C antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Integrin-linked kinase-associated serine/threonine phosphatase 2C protein (1-392AA)
Target	ILKAP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H0C8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.