



Product Datasheet

Keratin, type II cytoskeletal 4 antibody (orb241679)


Description	Keratin, type II cytoskeletal 4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Keratin, type II cytoskeletal 4 protein (1-534AA)
Target	KRT4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19013
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of A549 cells...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...